

The **multiSUB™ miniONE** electrophoresis system includes all of the equipment that you need to get up and running: gel tank, power supply and two casting sets. This gel tank/power supply combination is compact and easy to use.

multiSUB miniONE is an all in one horizontal electrophoresis unit. Featuring a built in power supply with voltage options of 35V, 50V and 100V, the miniONE is a versatile system suitable for a wide range of applications. A built in timer for runs from 0-99 minutes means you can set up your parameters and leave the system to complete the run automatically without fear of loosing bands.

The system comes complete with 2 gel casters for wide and mini gels as well as reversible combs for high throughput or high sample volume.

## KEY FEATURES

miniONE electrophoresis system is ideal for personal use, small laboratories or the classroom.

- All in one horizontal electrophoresis system
- In built power supply with 35V, 50V and 100V settings
- Timer function for runs from 0 99 minutes
- 2 gel sizes and reversible comb options





TECHNICAL SPECIFICATIONS		
Input Power	AC100~120V, 50~60Hz / AC200~240V, 50~60Hz	
Output Power	DC35V / DC50V / DC100V	
Gel Dimensions	Small: 54 x 60mm Large: 110 x 60mm	
Volume of Tank	230ml	
Construction of Bath	PC+ABS with high temperature resistance	
Timer range	1-99min	
Bath Dimensions	120 x 110 x 45mm	









pour and cast gel

place tray in unit and cover with buffer

load samples

start the run

Ordering Information			
MSMINIONE	Includes multiSUB™ miniONE electrophoresis system with	MSO-GCL	multiSUB miniONE Gel Caster – Large
	Built-in power supply,	MSO-GCS	multiSUB miniONE Gel Caster – Small
	2 x MSO-UVL, 4 x MSO-UVS, 1 x MSO-GCL, 1 x MSO-GCS,	MSO-UVL	multiSUB miniONE Large Gel Tray 110mm×60mm
	2 x MSO-1-5/9DS, 2 x MSO-1-12/22DS	MSO-UVS	multiSUB miniONE Small Gel Tray 54mm×60mm
MSO-1-12/22DS	Full Length Combs for miniONE Large Gel Tray		
MSO-1-5/9DS	Double Comb for miniONE Small Gel Tray		